

Specification Sheet

PRODUCT: Phosphate Buffered Saline-10x Concentration

| | | | |
|---------------------|------------------|------------------|-------------------------------------|
| CATALOG NO.: | 472842 or 492842 | LOT NO.: | TBD |
| MANUFACTURE DATE | TBD | EXPIRATION DATE: | 1 year from the date of manufacture |
| STORAGE CONDITIONS: | 4 - 40°C | | |

Test Results

| Test | Specification | Test Method | Result |
|-----------------------|------------------------|--------------------|--------|
| Appearance | Clear/Colorless Liquid | Visual | TBD |
| pH | 6.6 - 6.9 | Potentiometric | TBD |
| Heavy Metals(as Lead) | ≤ 5ppm | ICP | TBD |
| Cell Culture | Pass | Mouse Embryo Assay | TBD |
| Protease | None Detected | Enzymatic | TBD |
| DNase | None Detected | Fluorometric | TBD |
| RNase | None Detected | Fluorometric | TBD |
| Peroxides Content | ≤ 1ppb | Enzymatic | TBD |
| Sterility | No Growth | USP <71> | TBD |
| Endotoxins | ≤ 1.0 EU/mL | USP <85> | TBD |

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

Simple Lifecycle