

## Specification Sheet

**PRODUCT: Sodium Hydroxide, 1.0M**

CATALOG NO.:	472844 and 492844	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	1 Year from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C		

### Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
pH	≥ 13	Potentiometric	TBD
Identification (Sodium)	Pass/Fail	Flame Test	TBD
Protease	None Detected	Enzymatic	TBD
DNase	None Detected	Fluorometric	TBD
RNase	None Detected	Fluorometric	TBD
Peroxide	≤ 3ppb	Enzymatic	TBD
Assay	0.9 - 1.1N	Titration	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxins	≤ 2 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

### Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

*Simple Lifecycle*