

## Specification Sheet

**PRODUCT: Tris Buffered Saline-10X Concentration**

CATALOG NO.:	472845 and 492845	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	1 Year from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C		

**Test Results**

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
pH	7.4 - 7.8	Potentiometric	TBD
Osmolarity (1:10 Dilution)	280 ± 20 mOsm	Freeze Point Depression	TBD
Protease	None Detected	Enzymatic	TBD
DNase	None Detected	Fluorometric	TBD
RNase	None Detected	Fluorometric	TBD
Peroxide	≤ 1ppb	Enzymatic	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxins	≤ 1 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

**Manufacturing Site**

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538