

## Specification Sheet

**PRODUCT: Water for Reproduction**CATALOG NO.: 472847 and  
492847

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date  
of manufacture

STORAGE CONDITIONS: 4 - 40°C

**Test Results**

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
Total Organic Carbon	≤ 8.0 ppm	USP <643>	TBD
Conductivity	≤ 5.0 μS/cm @ 25°C	USP <645>	TBD
Protease	None Detected	Enzymatic	TBD
DNase	None Detected	Fluorometric	TBD
RNase	None Detected	Fluorometric	TBD
Peroxides Content	≤ 1ppb	Enzymatic	TBD
Cell Culture	Pass	Mouse Embryo Assay	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxins	≤ 0.25 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

**Manufacturing Site**

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

*Simple Lifecycle*