

## Specification Sheet

**PRODUCT: Hydrochloric Acid, 0.6N**

CATALOG NO.: 492446

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date  
of manufacture

STORAGE CONDITIONS: 4 - 40°C

PACKAGE SIZE: 3L

### Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
Chloride Ion	0.595 - 0.605N	Titration	TBD
Sterility	None Detected	USP <71>	TBD

Meets current USP/NF requirements.

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

### Manufacturing Site

Lifecycle Biotechnologies, 421 Marti Drive, Cleburne, TX 76033