

Specification Sheet

PRODUCT: Hydrogen Peroxide, 3%

CATALOG NO.: 492113

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date
of manufacture

STORAGE CONDITIONS: 4 - 40°C

PACKAGE SIZE: 4L

Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
pH	4.2 ± 0.5	Potentiometric	TBD
Assay	Direct Titration	3.0 ± 0.5%	TBD
Sterility	None Detected	USP <71>	TBD

Meets current USP/NF requirements.

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 421 Marti Drive, Cleburne, TX 76033