

Specification Sheet

PRODUCT: 0.9% Sodium Chloride (w/v)

CATALOG NO.: 492204

LOT NO.: TBD

MANUFACTURE DATE: TBD

EXPIRATION DATE: 1 Year from the date
of manufacture

STORAGE CONDITIONS: 4 - 40°C

PACKAGE SIZE: 4L

Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
Assay	0.85 - 0.95%	Titration	TBD
Endotoxins	≤ 0.5 EU/mL	USP <85>	TBD
Sterility	None Detected	USP <71>	TBD

Meets current USP/NF requirements.

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 421 Marti Drive, Cleburne, TX 76033