

Specification Sheet

PRODUCT: Water, Sterile Filtered

[Meets Testing Specification for USP Water for Irrigation]

CATALOG NO.:	492133	LOT NO.:	TBD
MANUFACTURE DATE	TBD	EXPIRATION DATE:	1 Year from the date of manufacture
STORAGE CONDITIONS:	4 - 40°C	PACKAGE SIZE:	4L

Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
Total Organic Carbon	≤ 8.0 ppm	USP <643>	TBD
Conductivity	≤ 5.0 μS/cm @ 25°C	USP <645>	TBD
Sterility	None Detected	USP <71>	TBD
Endotoxins	≤ 0.25 EU/mL	USP <85>	TBD

Meets current USP/NF requirements.

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 421 Marti Drive, Cleburne, TX 76033